



(<https://gem.gov.in/>)

[PRODUCTS](#) [SERVICES](#) [CONTENT](#)

Enter Keywords?



HOME ([HTTPS://MKP.GEM.GOV.IN/DASHBOARD](https://mkp.gem.gov.in/dashboard)) / BID FINALIZATION

BID DETAILS

TECHNICAL EVALUATION

FINANCIAL EVALUATION

EVALUATION

BID AWARDED

1. Bid Details

Your corrigendum has been published successfully. (<https://bidplus.gem.gov.in/buyer-bid-finalization/9407874>)

Bid Number: **GEM/2026/B/7610134** (/showbidDocument/9407874)

Bid Status: Active

Quantity: 7

Bid Validity (From End Date): 180 (Days)

Competent Authority Document:
View

Bid Start Date / Time: 19-06-2026 16:06:15

Bid End Date / Time: 18-09-2026 14:00:00

Bid Opening Date / Time: 19-09-2026 14:00:00

Consignees / Reporting Officer / Delivery Location(S)

EMD: Required Track EMD ([Https://Bidplus.Gem.Gov.In/Bidding/Track/Trac](https://bidplus.gem.gov.in/Bidding/Track/Trac))

Average Turn Over of Last 3 Years: 8 Lakh (s) **Experience with Gov. Required:** Year (s)

Buyer Details

Name: Poojan Marwaha

Address: Poojan Marwaha,Buyer77.Aiimsa.Hp@Gem buyer.In,Aiims ,Kothipura,Changar Palasia,Bilaspur,HIMACHAL PRADESH,174001,India

Ministry: Ministry Of Health And Family Welfare

Department: Department Of Health And Family Welfare

Organisation: All India Institute Of Medical Sciences (Aiims)

Office: Bilaspur Himachal Pradesh

Corrigendum Details

Modified On: 2026-09-08 12:44:23

[Download](#) (/bidding/buyer/showcorrigendumpdf/4567890/9407874)

Hash Value (Algorithm - SHA256):
94c6765c86ff384e6e7429b4d8f5988196ebd94555465a7d8ade06a5306d13c4

Modified On: 2026-09-08 12:38:01

Bid extended to **2026-09-18 14:00:00**

[Cancel Bid](#)

[Extend Bid](#)

[Edit Terms](#)

[Edit Pre Bid](#)

WEB INFO

[Terms of Use \(https://gem.gov.in/termsCondition\)](https://gem.gov.in/termsCondition)

[Website Policies \(https://gem.gov.in/websitePolicies\)](https://gem.gov.in/websitePolicies)

[Document Help \(https://gem.gov.in/help\)](https://gem.gov.in/help)

[Site Map \(https://gem.gov.in/sitemap\)](https://gem.gov.in/sitemap)

[Web Information Manager \(https://gem.gov.in/web-information-manager\)](https://gem.gov.in/web-information-manager)

ABOUT GeM

[Introduction to GeM \(https://gem.gov.in/aboutus\)](https://gem.gov.in/aboutus)

[Statistics \(https://gem.gov.in/statistics\)](https://gem.gov.in/statistics)

[Right to Information \(https://gem.gov.in/RTI\)](https://gem.gov.in/RTI)

[Analytics \(https://sso.gem.gov.in/ARXSSO/oauth/login\)](https://sso.gem.gov.in/ARXSSO/oauth/login)

[New on GeM \(https://gem.gov.in/new-categories\)](https://gem.gov.in/new-categories)

[BRAND GeM \(https://gem.gov.in/brand-gem\)](https://gem.gov.in/brand-gem)

NEWS & EVENTS

[Newsroom \(https://gem.gov.in/media\)](https://gem.gov.in/media)

[Gallery \(https://gem.gov.in/gallery\)](https://gem.gov.in/gallery)

[Notifications \(https://gem.gov.in/landing/index/allnews\)](https://gem.gov.in/landing/index/allnews)

[CCM Schedule \(https://gem.gov.in/ccm_data\)](https://gem.gov.in/ccm_data)

[Forums \(https://gem.gov.in/forum\)](https://gem.gov.in/forum)

TESTIMONIALS
([HTTPS://GEM.GOV.IN/TESTIMONIAL-BUYERS](https://gem.gov.in/testimonial-buyers))

RESOURCES

[GeM Handbook \(https://assets-bg.gem.gov.in/resources/pdf/GeM_handbook.pdf\)](https://assets-bg.gem.gov.in/resources/pdf/GeM_handbook.pdf)

[OM's/Circulars \(https://gem.gov.in/support/government_oms_circulars\)](https://gem.gov.in/support/government_oms_circulars)

[Terms and Conditions \(https://gem.gov.in/support/terms_conditions\)](https://gem.gov.in/support/terms_conditions)

[Policies/Manuals \(https://gem.gov.in/support/buyers\)](https://gem.gov.in/support/buyers)

[Miscellaneous \(https://gem.gov.in/support/miscellaneous\)](https://gem.gov.in/support/miscellaneous)

[MoUs \(https://gem.gov.in/mou\)](https://gem.gov.in/mou)

TRAINING

[LMS \(https://lms.gem.gov.in/\)](https://lms.gem.gov.in/)

[Training Calendar \(https://gem.gov.in/training\)](https://gem.gov.in/training)

[Training Module \(https://gem.gov.in/training/training_module\)](https://gem.gov.in/training/training_module)

[Facilitators \(https://gem.gov.in/training/facilitators\)](https://gem.gov.in/training/facilitators)

[Download GeM Logo \(https://gem.gov.in/resources/images/download-gem-logo.jpg\)](https://gem.gov.in/resources/images/download-gem-logo.jpg)

NEED HELP?

[FAQs \(https://gem.gov.in/userFaqs\)](https://gem.gov.in/userFaqs)

[Video Guides \(https://gem.gov.in/training/videos/buyers\)](https://gem.gov.in/training/videos/buyers)

[Raise-a Ticket \(https://gem.gov.in/gemtickets/create\)](https://gem.gov.in/gemtickets/create)

[Contact Us \(https://gem.gov.in/contactUs\)](https://gem.gov.in/contactUs)

[Careers \(https://gem.gov.in/landing/index/careers\)](https://gem.gov.in/landing/index/careers)





Revised Technical Specifications
Multipara Monitor with Central Monitoring System (CMS) – 8th August 2026



1. Central stations must display comprehensive patient information from bedside monitors in a simple and effective manner for in-depth analysis and to increase workflow productivity.
2. Maximum number per display: 8 on a single display.
3. Maximum number of bedside monitors/Central Station: ≥ 16 .
4. Central Display: LCD.
5. Type of display: Single.
6. Display size: Minimum 19 inches.
7. Maximum number of main display (patient sector) selectable: 16.
8. Maximum number of individual patient display: 2.
9. Trending: Graphical, minimum 72-hour trend storage.
10. ST Analysis: Yes.
11. Event storage: Minimum 120 events.
12. Alarm default settings for low/high BP, low saturation, low/high respiration with volume setting.
13. Hardware platform: Closed system.
14. Operating system: Windows.
15. System security: Yes.
16. Maximum number of devices on a single network: ≥ 16 .
17. Bedside monitor-to-bedside monitor communication.
18. Central ECG data management system.
19. Strip chart recorder: 2 channels.
20. Printer: Laser printer with UPS (30-minute backup).



9/9/26



21. Warranty: 5 years.
22. User/Technical/Maintenance manuals in English (hard and soft copy).
23. Demonstration and training after installation before acceptance.
24. Principal Manufacturer must have direct presence/approved service centre in India.
25. Copies of reports/certifications to be furnished on demand.
26. Attach certificate from Principal Manufacturer confirming CMS software support for next 10 years with free updates if required.
27. Type of Monitor: Preconfigured/Non-modular.
28. Display: Minimum 12.1-inch colour LCD touchscreen.
29. Resolution: 800×600 pixels or better.
30. Rechargeable battery with more than 60 minutes operating time.
31. 12-lead ECG interpretation.
32. Adult Heart Rate (ECG): 10–300 beats/min.
33. ST segment range: -2.0 mV to $+2.0$ mV.
34. Pulse Rate (SpO₂): 30–240 beats/min with ± 3 bpm accuracy.
35. BP Range: Adult 10–270 mmHg; Child 10–235 mmHg; Neonate 10–135 mmHg; Mean 20–235 mmHg.
36. Repeat of Point 34.
37. Cuff pressure range: 0–280 mmHg.
38. Probe type: Reusable.
39. Cuffs supplied: Adult (Medium), Paediatric (Standard), Neonatal (Standard).
40. Temperature probe type: Skin.
41. Temperature range: 0°C – 45°C or more with 80–100% accuracy.
42. SpO₂ range: 1–100% with 80–100% accuracy.
43. Auto measurement interval: 1,2,3,4,5,10,15,20,30,45 minutes.
44. Duration of alarm silence: 120 seconds.



45. Connectivity provision with Central Nursing Station.
46. Reusable SpO₂ probes: Adult, Child, Neonate.
47. Accessories for wall mounting.
48. Battery type: Lead Acid / Nickel-metal hydride / Lithium battery.
49. Battery backup time: 60 minutes.
50. Battery recharge time: ≤ 6 hours.
51. Toll-free technical support by OEM/authorised agency.
52. Warranty: 5 years.
53. Accessories: ECG leads-5; ECG cable-1; Disposable electrodes-1 packet; NIBP cuffs (Adult, Paediatric, Neonatal)-3 sets each; Temperature probes (Adult, Paediatric, Neonatal)-1 each; SpO₂ probes (Adult, Paediatric, Neonatal)-3 each; Power cord-1; Rechargeable battery-1.

